

INDION® MSR

Mercury Specific Resin

Description

INDION MSR is a crosslinked macroporous ion exchanger having thiol as functional groups. It has been specifically developed for the selective adsorption of mercury and noble metals. INDION MSR is completely insoluble in water and all common solvents. It is stable in acidic as well as alkaline pH range. INDION MSR has excellent physical stability due to its macroporous structure. It is unaffected by reducing agents. However, INDION MSR can get decrosslinked in the presence of oxidising agents such as chlorine, peroxides, permanganates etc. It is therefore essential to remove all the oxidants from the process stream, before it comes in contact with INDION MSR.

Applications

INDION MSR is specifically desired for removing mercury from the waste water, where mercury exists in either ionic or elemental forms. Mercury can be brought down to 0.005 ppm from 20 ppm using INDION MSR. INDION MSR can also be used as a polisher for removing the last traces of mercury,

where a major quantity has already been removed by earlier processes. INDION MSR's efficiency in removing mercury in monovalent or elemental form is relatively low. Hence, oxidising mercury to the divalent state may be necessary. It is essential that any oxidising agent in this process be removed before passing through the resin column. INDION MSR should be preferably used in the acidic or neutral pH range. INDION MSR's capacity for removal of mercury is influenced by many factors. These include,

- Concentration of mercury in the feed stream
- Whether mercury exists as organic or inorganic compounds
- The valency of mercury
- pH of the solution
- Temperature of influent
- Permitted leakage of mercury in the effluent
- Service flow rate
- Regeneration level

INDION MSR's mercury removal capacity is not affected by common constituents of natural water such as hardness, alkalinity or inorganic salts.

However, mercury removal from cyanide solutions is not recommended.

Characteristics	
Appearance	Opaque pale yellow to brown beads
Matrix	Crosslinked polystyrene
Functional Group	Thiol
lonic form as supplied	Hydrogen
Weight capacity in H form	3.6 meq/dry g, minimum
Moisture holding capacity	38 - 43%
Shipping weight*	690 - 730 kg/m³
Particle size range	0.3 to 1.2 mm
> 1.2 mm	5.0%, maximum

< 0.3 mm	1.0%, maximum
Uniformity co-efficient	1.7, maximum
Effective size	0.40 to 0.55 mm
Maximum operating temperature	60° C
Operating pH range	1 to 10
Resistance to reducing agents	Good
Resistance to oxidizing agents	Generally good, chlorine should be absent
*Weight of resin, as supplied, occupying 1 m³ in a unit after backwashing and draining.	

Elution of mercury

30% HCl has been proved to be an efficient regenerant to elute the mercury adsorbed on INDION MSR. This helps the chlor-alkali industry to recycle the expensive mercury to the cell house thereby totally eliminating the problem of mercury disposal.

Packing

Fiber drums with liner bags

HDPE Lined bags : 25/50 lts
LDPE bags : 1 cft/25 lts
Super sack : 1000 lts
Super sack : 35/40/42 cft
MS/HDPE drums with liner bags : 180/200 lts

Storage

Ion exchange resins require proper care at all times. The resin must never be allowed to become dry. Regularly open the plastic bags and check the condition of the resin when in storage. If not moist, add enough clean demineralised water and keep it in completely moist condition. Always keep the resin drum in the shade.

Safety

Acid and alkali solutions used for regeneration are corrosive and should be handled in a manner that will prevent eye and skin contact. If any oxidising agents are used, necessary safety precautions should be observed to avoid accidents and damage to the resin.

INDION range of lon Exchange resins are produced in a state-of-the-art ISO 9001 and ISO 14001 certified manufacturing facilities at Ankleshwar, in the state of Gujarat in India.

7 cft

To the best of our knowledge the information contained in this publication is accurate. Ion Exchange (India) Ltd. maintains a policy of continuous development and reserves the right to amend the information given herein without notice.

INDION® is the registered trademark of Ion Exchange (India) Ltd.



ION EXCHANGE (INDIA) LTD.

Corporate Office

Ion House, Dr. E. Moses Road, Mahalaxmi, Mumbai - 400011 | Tel: +91 22 6231 2000 E-mail: ieil@ionexchange.co.in

Regional and Branch Offices

Bengaluru | Bhubaneswar | Chandigarh | Chennai Delhi | Hyderabad | Kolkata | Lucknow | Vadodara Vashi | Visakhapatnam

International Division

R-14, T.T.C MIDC, Thane - Belapur Road, Rabale, Navi Mumbai - 400 701 | Tel: +91 22 6857 2400 E-mail: export.sales@ionexchange.co.in

Overseas Offices

Bangladesh | Canada | Indonesia | Kenya Malaysia | Oman | Portugal | Saudi Arabia | Singapore South Africa | Sri Lanka | Tanzania | Thailand | UAE | USA

Manufacturing Units

India - Ankleshwar | Hosur | Patancheru | Rabale | Verna | Wada Overseas - Bangladesh | Indonesia | Saudi Arabia | UAE

All India Service and Dealer Network

